


photo logging and event insights

You can use the **Dexcom G7 app** to log daily events—like insulin doses, meals, activity, meter values, medications, and even quick notes. Tracking these events can help you understand how your daily actions affect your glucose.

In your **Dexcom G7 app**, events can be added in the **Log tab** or by tapping the **plus (+) button**.

PHOTO LOGGING

Make your event more visual by adding a photo when you log it. Add a picture of your meal, the exercise you did, medication, or any other type of event. This helps provide additional context when looking back at a particular event.



SMART FOOD LOG



Take a photo of food or upload it from the camera roll



AI identifies ingredients and populates meal description



Add notes like carbs, date and time

EVENT INSIGHTS

Events Insights help take event logging to the next level. In the **Log tab**, you can review logged events and see the glucose impact with a before/after trend graph and time stamps to better understand the “why” behind any glucose peaks and valleys.

For example:

- If you log a meal, you’ll see a snapshot of how your glucose changed from mealtime to three hours afterward, allowing you to see what happened in your body over that period.
- If you log an activity, you’ll see a snapshot of what your glucose did from start to end – for example, how it shifted during a 60-minute run.

