

## Comparing Continuous Glucose Monitoring (CGM) and Blood Glucose Monitoring (BGM) in Adults With Inadequately Controlled, Insulin-Treated Type 2 Diabetes: A 12-Month, Single-Center, Randomized Controlled Trial<sup>1</sup>

### PURPOSE OF STUDY

To compare the 12-month effects of **Dexcom RT-CGM** versus **BGM** in **adults with poorly controlled (A1C  $\geq 7.5\%$ ), insulin-treated type 2 diabetes.**

### DEMOGRAPHICS



Adults with T2D for over 1 year  
(n=76)  
( $\geq 18$  years of age)



HbA1c  $\geq 7.5\%$



Insulin at least once daily


**(56% were on at least 3 medications)**



Dexcom G6\*

### METHODS

12-month single-center, open-label, parallel, RCT



One of the first T2 RCTs using CGM that was conducted in a study population with a high use of GLP1RA and SGLT2i.



Participants in both groups received a 3-h diabetes self-management education (DSME).



Minimum study contacts replicated real-world usual care. HCP visits took place every 3 months.

\* CGM – Continuous glucose monitoring, BGM – Blood glucose monitoring, TIR – Time in range (70–180 mg/dL), TAR– Time above range (>250 mg/dL) The mean active sensor time assessed during the last blinded CGM at 12-month follow-up was 96.3%

† Group results from baseline to 12 months. ‡ (95% CI: 4.6; 25.9)% (P = 0.006) § (P = 0.006)

1. Lind N, et al. Diabetes Care. 2024;47(5):881–9 2. Battelino T, Danne T, Bergenstal RM, et al. Diabetes Care 2019;42(8):1593–1603; doi: 10.2337/dci19-0028

## GLYCEMIC AND METABOLIC OUTCOMES

When compared to BGM, CGM\* use was associated with significantly improved outcomes.<sup>†</sup> The improvements in A1C together with a decrease in total daily dose (TDD) of insulin, BMI and body weight was achieved without increasing other glucose-lowering treatments<sup>1</sup>

↑15.2%<sup>‡</sup>

MORE TIR<sup>2</sup>  
(3h 39min a day)

↓15.8%<sup>§</sup>




DECREASE IN TAR\*<sup>‡</sup>

↓10.6

UNITS LESS  
OF INSULIN

↓3.3 kg

DECREASE IN  
BODY WEIGHT  
(SELF-REPORTED)



**7.5%** Participants in the CGM group no longer required insulin for treatment at 12 months (p=0.13).<sup>1</sup>

## Key Takeaways for Dexcom CGM

- **Significant improvement of TIR<sup>2</sup> in the Dexcom RT-CGM treated group compared to BGM treatment at 6 months which sustained until 12 months.** (6 months: 12.4% TIR improvement (2 h 59 mins/day), P=0.021. 12 months: 15.2% TIR improvement (3 h 39 mins/day), P=0.006).<sup>1</sup>
- **The empowerment of people to make informed changes in their lifestyle** is supported by the between-group difference in TIR peaking during the day and changes in diet and activity.<sup>1</sup>
- **Treatment de-escalation is supported by significant reduction of total daily dose of insulin (TDD),** and a tendency to lower MDI and lower insulin dependency, without increase in other glucose-lowering treatments.
- **Dexcom RT-CGM significantly improved glycemic, and self-reported metabolic participant outcomes compared to BGM in adults with insulin-treated type 2 diabetes.**



**Minimal education was needed** to empower self-managed behavior change.<sup>1</sup>



The increase of TIR was achieved **without an increase in TBR.**<sup>1</sup>



Participants **self-reported greater improvements** in:

- General well-being
- Diabetes related distress
- Treatment and monitoring satisfaction.<sup>1</sup>

\* CGM – Continuous glucose monitoring, BGM – Blood glucose monitoring, TIR – Time in range (70-180 mg/dL), TAR – Time above range (>250 mg/dL) The mean active sensor time assessed during the last blinded CGM at 12-month follow-up was 96.3%  
<sup>†</sup> Group results from baseline to 12 months. ‡ (95% CI: 4.6; 25.9%) (P = 0.006) § (P = 0.006)  
<sup>1</sup> Nanna Lind, Merete B. Christensen, Dorte L. Hansen, and Kirsten Nørgaard, Diabetes Care 2024;47(5):1-9 | https://doi.org/10.2337/dc23-2194 2. Battelino T, Danne T, Bergenstal RM, et al. Diabetes Care 2019;42(8):1593-1603; doi: 10.2337/dci19-0028  
 Dexcom is a registered trademark of Dexcom, Inc. in the United States and/or other countries. ©2024 Dexcom, Inc. All rights reserved. MAT-3802 v1.0

**dexcom**