

Real-World Use of Control-IQ Technology is Associated with a Lower Rate of Severe Hypoglycemia and Diabetic Ketoacidosis Than Historical Data: **Results of the Control-IQ Observational (CLIO) Prospective Study¹**

Background

Individuals on intensive insulin therapy (IIT) face an elevated risk of experiencing **Severe Hypoglycemia (SH)** and **Diabetic Ketoacidosis (DKA)**. It is important to minimize the occurrence of these **adverse events (AE)**, when possible. Automated Insulin Delivery (AID) systems may help reduce some of these events.

Purpose of Study



To assess the safety and effectiveness of the Tandem t:slim X2 insulin pump with Control-IQ technology in the real-world outpatient setting over 12 months of use.

Participant Demographics (n=3157)

Single-arm, prospective, fully virtual, longitudinal study



T1D



≥ 6 y.o.



Control IQ & Dexcom
G6 for 12 months



Virtual check-in and
data uploads


Study Outcomes Measured

- ✓ The primary outcomes measured were incidence rates of:
 - SH and DKA
 - Safety of the automatic population of CGM readings into the bolus calculator when Control-IQ technology was in use.
- ✓ Glycemic outcomes as a measure of efficacy of the system.
- ✓ Patient-reported satisfaction with and trust in the system, diabetes impact, and sleep quality of Control-IQ with Dexcom CGM technology users.


Main Outcomes



Severe Hypoglycemia (SH)*



Diabetic Ketoacidosis (DKA)*



SH and DKA rates were lower for users of **t:slim X2 with Control-IQ with Dexcom CGM** compared to historical data[§] for both adults and children.



Average time in range was **>70%** for adults and **>60%** for children, both groups achieving international consensus guidelines¹



42% reported a decrease of **diabetes impact** on quality of life[†]

25%

increase in device satisfaction[‡] in both participants transitioning to Control-IQ from a different pump system and for participants transitioning from MDI therapy

Key Takeaways for t:slim X2 with Control IQ and Dexcom CGM

- **Lower observed rates of AE's** occurred independent of baseline A1C or prior insulin delivery method²
- Control-IQ with Dexcom CGM technology is associated with a **lower rate of severe hypoglycemia** (compared to historical data)²
- The study also indicates a **lower rate of diabetic ketoacidosis** in individuals using Control-IQ with Dexcom CGM technology compared to historical data²
- **Reduction in diabetes burden** was consistently reported²

*Mean (SD) AE rates per 100 patient years for study participants as self-reported on monthly surveys through 12 months of Control-IQ technology use, compared to historical rates

AE = Adverse events

† Diabetes Impact median score was 4.75 (3.25–6.0) at baseline and 2.75 (2.0–4.0) at 12 months (P < 0.01).

‡ DIDS Device Satisfaction median score increased from 7.29 (5.57–8.71) at baseline to 9.14 (8.29–9.71) at 12 months

§ Data from T1D Exchange Registry

1. Battelino T, Danne T, Bergenstal RM, et al. Clinical targets for continuous glucose monitoring data interpretation: Recommendations from the International Consensus on Time in Range. *Diabetes Care* 2019;42(8):1593–1603; doi: 10.2337/dci19-0028

2. Graham R, et al. *Diabetes Technol Ther*. 2024 Jan;26(1):24–32. doi: 10.1089/dia.2023.0341.

For more information, visit provider.dexcom.com

Dexcom is a registered trademark of Dexcom, Inc. in the United States and/or other countries. ©2024 Dexcom, Inc. All rights reserved. MAT-2235 v1.0